

In Vivo Transfection

In Vivo Transfection Reagents and Delivery Protocols

IN VIVO [PEG LIPOSOME](#) NANOPARTICLE LIPID-BASED POLYMER-BASED USES OF DNA
TRANSFECTION [TRANSFECTION](#) TRANSFECTION TRANSFECTION TRANSFECTION TRANSFECTION

PEG LIPOSOME TRANSFECTION

PEGYLATED LIPOSOMES

In order to enable intracellular delivery of drugs and medicines, a vehicle molecule is required. Such vehicle molecules vary in type, mechanism of action, targeting of specific organs or tissues, and unique physical and chemical properties. PEG-liposomal delivery technology is known to be superior in terms of delivery efficiency and tissue targeting applications while inducing minimal immune response. Preclinical contract research organizations (such as [Altogen Labs](#)) provide PEG liposome encapsulation services (see service description [here](#)) for *in vivo* [RNAi](#) and multiple other gene expression *in vivo* applications. Encapsulated mRNA, siRNA, shRNA, plasmid DNA, proteins and small molecules are used for both *in vitro* and *in vivo* transfection laboratory experiments.

Cationic polymers that form stable positively charged nanoparticles with DNA in water are commercially available for [in vivo transfection](#). Cationic lipids are commonly used to deliver nucleic acids due to electrostatic interactions with amphiphilic molecules. There are many methods of introducing nucleic acids into an organism usually referred to as [gene therapy](#). Some of the more successful ones include:

- • Injection of naked DNA – mixed success, often cleared without effect
- • Sonoporation [method](#)
- • Electroporation [method](#)
- • Lipoplexes [method](#)
- • Polyethylenimine (PEI) reagent can be used to deliver genes or RNA via a spinal catheter *in vivo* for successful [transfection](#).

PEG (POLYETHYLENE GLYCOL) TRANSFECTION COMPLEXES


A liposome is a very simple cellular structure that can carry other molecules inside the cell, including small drugs or a segment of genetic code (DNA or RNA). These fundamental structures are found abundantly in and around cells and can be artificially created. Historically, liposome mediated delivery has been utilized for many decades. However, *in vivo* liposome delivery systems experience immediate elimination from circulation. Positively charged complexes interact with negatively charged components in plasma, with liposome binding to lipoproteins resulting in inactivating the complexes. The lipoprotein interaction with lipoplexes forms a bond between the two components that leads to aggregate formation.

To circumvent elimination and aggregation, PEG-lipids can be inserted into the bilayer of the lipid complex. The addition of the PEG-lipids shields the lipid complex from being cleared from circulation to increase bioavailability and limits aggregation.

PEG-lipids are commonly used in pharmaceuticals due to its ability to stabilize liposomes and low molecule cost. Formulation chemists routinely alter the PEG chain length, PEG density, spacers and linkers to customize zeta-potential, complex size and plasma interaction.

PEG-liposomes are liposomal complexes (cargo molecules include plasmid DNA, mRNA, small RNA, small molecules, and proteins) that do not induce immune response and provides superior *in vivo* delivery compared to many

counterparts. PEG-liposomes are commercially available from Altogen Biosystems – [PEG-Liposome Kits](#).



PEG IN VIVO TRANSFECTION PROTOCOL PEGYLATED LIPOSOME DESCRIPTION

PEGylated cationic lipid liposome-based *In Vivo* Transfection Reagent is a lipid-functionalized PEG complex proprietary formulation optimized for *in vivo* delivery of miRNA, siRNA, plasmid DNA and proteins. PEG-Liposome complexes are stable for at least 16 hours in serum. There is no detectable inflammatory response or toxicity due to the PEG modification. Efficient delivery has been exhibited to these tissues: pancreas, spleen, kidney, liver and tumors via systemic administration.

PEGYLATED LIPOSOME SYSTEMIC ADMINISTRATION (I.V. INJECTION)

Below is a recommended protocol for the intravenous injection of PEGylated liposomes for *in vivo* experiments.

1. Dilute 60 μg of plasmid DNA or 100 μg of siRNA in 100 μL nuclease-free water and vortex gently.
2. Add 100 μL of the diluted to a sterile tube containing 50 μL Transfection Reagent.
3. Incubate for 15-20 min at room temperature.
4. Add 10 μL of Transfection Enhancer Reagent and vortex gently to mix.
5. Incubate for 5 min at room temperature.
6. Add required amount of sterile solution of 5% glucose (w/v):

Animal body weight (g)	Final injection volume (ml)
10-14	0.2
15-19	0.3
20-24	0.4
25-29	0.5
30-35	0.6

7. Inject animals: Please note that delivery efficiency can be significantly increased by performing a secondary injection (at least 12 hours after first injection).
8. Maximum mRNA target expression effect is typically observed 12-36 hours after injection, while maximum effect on protein level is achieved 24-48 hours post-injection.




Figure: Systemic administration (i.v. injection) of PEGylated liposome In Vivo reagent conjugated with siRNA targeting Lamin A/C mRNA or non-silencing negative control siRNA. Tissues were collected and RNA was isolated 48 hours after first injection. Samples were analyzed by qPCR for Lamin A/C gene expression levels and normalized to a housekeeping gene. Data are means \pm SD (n=6).

